

Lyofast MW 039 R

Description

Lyofast MW 039 R consists of undefined strains of *Lactococcus lactis* ssp. *lactis*, *Lactococcus lactis* ssp. *cremoris* and selected strains of *Lactococcus lactis* ssp. *lactis* biovar *diacetylactis*, and *Leuconostoc* ssp.. Lyofast MW 039 R ensures a uniform and controlled production of sour cream like Smetana, fermented milk, fresh cheese, soft cheese, and semi-hard cheese. Lyofast MW 039 R has fast citrate fermentation.

Application

Sprinkle the culture powder directly into process milk under aseptic conditions ensuring that the culture is well dispersed by gentle stirring. The following may be used as inoculation guidelines:

Product	UC/100 l	Product	UC/100 l
Fresh cheese	0.5-2.0	Soft cheese	0.7-2.0
Semi-hard cheese	1.0-4.0	Fermented milk	0.5-2.0
Sour cream/Crème fraîche	1.0-4.0		

Rotation


The recommended rotation available upon request.

Acidification information

Standardised laboratory acidification test is conducted in milk powder, reconstituted at 9%, at defined temperature.

Acidification profile: inoculation level corresponding to 1 UC per 100 litres milk.

Standard activity: expressed as temperature/time/pH relations: 32°C/6.5 hours/pH 5.2 ± 0.1.



Culture information

Data are obtained under standardised laboratory conditions, and consequently, should be considered as guidelines.

Optimal temperature for growth	25-35 °C	Scalding temperature	Max. 39°C
Acidification capability	pH 4.5	Gas production/citrate/urea	++
Diacetyl production	+++		

Lyofast MW 039 R

Storage	Unopened pouches should be kept below -17°C.		
Package data	The freeze-dried culture is packed in waterproof and airproof aluminium pouches. The packaging material is food grade. Lyofast MW 039 R is available in 10 and 50 UC.		
Shelf life	18 months when stored at or below -17°C.		
Heavy metal specification	Pb (lead)	< 1 ppm	
	Hg (mercury)	< 0.03 ppm	
	Cd (cadmium)	< 0.1 ppm	
	<i>* Level of heavy metals is controlled on regular basis.</i>		
Microbiological specification	<i>Bacillus cereus</i>	<100 CFU/g	Method: Sacco M10 (1)
	Coagulase positive staphylococci*	<10 CFU/g	Method: Sacco M11(2)
	Enterobacteriaceae	<10 CFU/g	Method: Sacco M2 (3)
	<i>Escherichia coli</i>	<1 CFU/g	Method: Sacco M27 (4)
	<i>Listeria monocytogenes</i> *	Not detected in 25 g	Method: Sacco M13 (5)
	Moulds & yeasts	<10 CFU/g	Method: Sacco M3 (6)
	<i>Salmonella spp</i> *	Not detected in 25 g	Method: Sacco M12 (7)
	<i>* Analysed on regular basis. All analytical methods are available upon request.</i>		
	<i>(1)ISO 7932; (2)ISO 6888-1-2; (3)ISO 21528-1-2; (4)ISO11866-1-2/IDF 170-1-2; (5)ISO 11290-1-2; (6)ISO 6611/IDF 94; (7)ISO 6785/IDF 93.</i>		
GMO	The microbial strains are not genetically modified (GMO) in accordance with the European Directive 2001/18/EC. The strains are isolated from natural sources. The raw materials used are also GMO free in accordance with Regulation (EC) No. 1829/2003 and Regulation (EC) No.1830/2003. Statement available upon request.		
Allergens	The raw materials used are generally based on dairy ingredients. All materials are free of the following components and their derivatives: peanut, tree nut, sesame, egg, fish, shellfish, mollusc, crustacean, sulphite, cereals containing gluten, celery, mustard, soy and lupine. Statement available upon request.		
Safety information	Material Safety Data Sheet available on www.saccosrl.it		
Certificate	Lot certificate available upon request.		
ISO Kosher approval	Sacco S.r.l. is UNI EN ISO 9001:2008 certified since 1998, ISO 22000:2005 and FSSC 22000 certified since 2014. Sacco cultures are generally Kosher approved except for surface ripening cultures.		
Service	Please contact your distributor for guidance and instructions for your choice of culture and processing. Information about additional package sizes and sales units is also available upon request.		
Liability	This information is based on our knowledge trustworthy and presented in good faith. No guarantee against patent infringement is implied or inferred.		